

B. Supplementary Tables

Supplementary Table 1: List of antibodies

Antibody	Cat. No	Species	Application	Dilution
GFP	A10262 (Invitrogen)	Chicken	IHC-Fr	1:300
HNF4 α	sc-6556 (Santa Cruz)	Goat	IHC-Fr	1:50
CK19	602-670 (Abdomax)	Rabbit	IHC-Fr	1:100
CD31	557355 (BD Pharmingen)	Rat	IHC-Fr	1:100
CD45	NB100-77417 (Novus)	Rat	IHC-Fr	1:100
α -SMA	19245S (Cell Signaling)	Rabbit	IHC-Fr	1:300
F4/80	sc25830 (Santa Cruz)	Rabbit	IHC-Fr	1:100
HGF	AF-294-NA (R&D)	Goat	IHC-Fr	1:20
EGF	MAB236-100	Mouse	IHC-Fr	1:30
CD31-PE/Dazzle 594	102429 (BioLegend)	Rat	FC	1:40
CD11b-PE efluor 610	61-0112-80 (Invitrogen)	Rat	FC	1:80
CD45-PE efluor 610	61-0451-80 (Invitrogen)	Rat	FC	1:40

IHC-Fr: Immunohistochemistry frozen fixed tissues; FC: Flow Cytometry

Supplementary Table 2: Sequences of nucleoside-modified mRNA

Luciferase	ATGGAGGACGCCAAGAACATCAAGAAGGGCCCCGCCCTTCTACCCCTGGAGGACGGCACC GCCGGCGAGCAGCTGCACAAGGCC ATGAAGCGgTACGCCCTGGTGCCCGGCACCATCGCCTTACCCGACGCCACATCGAGGTGGACATCACCTACGCCGAGTACTTCGAGAT GTCCGTGCGCCTGGCCGAGGCCATGAAGCGgTACGGCCTGAACACCAACCACCGCATCGTGGTGTGCTCCGAGAACTCCCTGCAGTTCT TCATGCCCGTGTGGGCGCCTGTTTCATCGGCGTGCCGTGGCCCCGCCAACGACATCTACAACGAGCGGAGCTGCTGAACCTCCATG GGCATCTCCAGCCACCGTGTTGTTCTGTGTCGAAGAAGGGCCTGCAGAAAGATCCTGAACGTGCAGAAAGAAGCTGCCCATCATCCAGA AGATCATCATCATGGACTCCAAGACCGACTACCAGGGCTTCCAGTCCATGTACACCTTCGTGACCTCCACCTGCCCGCGGCTTCAACG AGTACGACTTCGTGCCGAGTCTTCGACCGCGACAAGACCATCGCCCTGATCATGAACCTCCTCGGCTCCACCGGCTGCCAAGGGC GTGGCCCTGCCACCGCACCCTGCGTGCCTTCTCCACGCCGCGACCCCATCTTCGGAACAGATCATCCCCGACACCGCCATC CTGTCCGTGGTGCCTTCCACACGGCTTCGGCATGTTACCACCTGGGCTACCTGATCTGCGGCTTCGCGTGGTGTGATGTACCGC TTCGAGGAGGAGCTGTTCTGCGCTCCTGCAGGACTACAAGATCCAGTCCGCCCTGCTGGTGCCACCTGTTCTCTTCTTCGCCAAG TCCACCTGATCGACAAGTACGACCTGTCCAACCTGCACGAGATCGCCTCCGGCGGCGCCCCCTGTCCAAGGAGGTGGGCGAGGCCG TGGCCAAGCGgTTCACCTGCCCGCATCCGCCAGGGCTACGGCTGACCGAGACCACCTCCGCCATCCTGATCACCCCGAGGGCGAC GACAAGCCCGCGCGCTGGGCAAGGTGGTGCCCTTCTCGAGGCCAAGGTGGTGGACCTGGACACCGGCAAGACCTGGGCGTGAAC CAGCGCGCGAGCTGTGCGTGCGCGGCCCATGATCATGTCCGGCTACGTGAACAACCCGAGGCCACCAACGCCCTGATCGACAAGG ACGGCTGGTGTCACTCCGGCGACATCGCTACTGGGACGAGGACGAGCACTTCTTCATCGTGACCGCTGAAGTCCCTGATCAAGTAC AAGGGCTACCAGGTGGCCCCCGGAGCTGGAGTCCATCTGCTGCAGCACCCCAACATCTTCGACGCGCGCTGGCCGGCTGCCG ACGACGACGCCGCGAGCTGCCCGCGCGTGGTGGTGTGGAGCACGGCAAGACCATGACCGAGAAGGAGATCGTGGACTACGTGG CCTCCAGGTGACCACCGCCAAGAAGCTGCGCGCGCGTGGTGTTCGTGGACGAGGTGCCCAAGGGCCTGACCGGCAAGCTGGACG CCCGCAAGATCCGCGAGATCCTGATCAAGGCCAAGAAGGGCGGCAAGATCGCCGTG
eGFP	ATGGTGAGCAAGGGCGAGGAGCTGTTACCGGGGTGGTGCCATCCTGGTTCGAGCTGGACGGCGACGTAAACGGCCACAAGTTCAGC GTGTCCGCGGAGGGCGAGGGCGATGCCACTACGGCAAGCTGACCTGAAGTTCATCTGCACCACCGGCAAGCTGCCGTGCCCTGGC CCACCCTCGTGACCACCTGACCTACGGCGTGCACTGCTTACGCGCTACCCCGACCATGAAGCAGCAGACTTCTCAAGTCGGCCA TGCCCGAAGGCTACGTCCAGGAGCGACCATCTTCTTCAAGGACGACGGCAACTACAAGACCCGCGCGAGGTGAAGTTCGAGGGCGA CACCTGGTGAACCGCATCGAGCTGAAGGGCATCGACTTCAAGGAGGACGGCAACATCCTGGGGCACAAGCTGGAGTACAAC AGCCACAACGTCTATATCATGGCCGACAAGCAGAAGAAGCGCATCAAGGTGAACCTTCAAGATCCGCCACAACATCGAGGACGGCAGCG TGCAGCTCGCCGACCACTACCAGCAGAACACCCCATCGGCGACGGCCCCGTGCTGCTGCCGACAACCACTACCTGAGCACCAGTCC GCCCTGAGCAAAGACCCCAACGAGAAGCGGATCACATGGTCTGCTGGAGTTCGTGACCGCCGCGGGATCACTCTCGGCATGGACG AGCTGTACAAG
HGF	ATGTGGGTGACCAAGCTGTGCCCCCTGCTGCTGCAGCAGTGTGCTGACCTGTGCTGCTGCCATCGCCATCCCTACGCCGA GGGCCAGCGCAAGCGCCGCAACACCATCCACGAGTTCAGAAGTCCGCAAGACACCTGATCAAGATCGACCCCGCCTGAAGATC AAGACCAAGAAGGTGAACACCGCCGACAGTGCGCCAACCGCTGCACCGCAACAAGGGCCTGCCCTTCACTGCAAGGCCTTCGTGTT CGACAAGGCCGCAAGCAGTGCCTGTGGTTCCTTCACTCCATGTCTCCGCGGTGAAGAAGGAGTTCGGCCACGAGTTCGACCTGT

[illegible]